GenePharma

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吉玛基因
www.genepharma.com

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ARTICLE

Open Access

Antisense oligonucleotides targeting lncRNA AC104041.1 induces antitumor activity through Wnt2B/ β -catenin pathway in head and neck squamous cell carcinomas

Mengwei Li^{1,2}, Xu Ding³, Yinan Zhang^{1,2}, Xin Li^{1,2}, Haoze Zhou^{1,2}, Li Yang^{1,2}, Yilin Li^{1,2}, Peiwei Yang^{1,2}, Xiaomin Zhang³, Jialiang Hu^{1,2}, Edouard Nice⁴, Heming Wu³ and Hanmei Xu^{1,2}

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
GenePharma

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www.genepharma.com

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- a. Cell proliferation assay of SCC4 cells transfected with AC104041.1 specific LNA-ASO or combined with salinomycin treatment (n = 5). Data are presented as the mean values ± SEM, the experiment was repeated three times, ***P < 0.001, compared with control cells (two-way ANOVA).
- b. Cell migration assay of SCC4 cells transfected with AC104041.1 specific LNA-ASO or combined with salinomycin treatment. Data are presented as the mean values ± SD, the experiment was performed in triplicates and repeated three times, **P < 0.01, ***P < 0.001 (Student's t-test).
- c. Tumor growth in Balb/c nude mice and representative immunohistochemical images of Ki-67 from tumors after subcutaneous injection of AC104041.1 specific LNA-ASO or combined with salinomycin treatment (n = 5 mice for each group). Data are presented as the mean values ± SEM, **P < 0.01, ***P < 0.001 (two-way ANOVA).
- d. Immunofluorescence of β-catenin (green) and nuclei (blue) in SCC4 cells transfected with AC104041.1 specific LNA-ASO or combined with salinomycin treatment. Scale bar represents 20 μm.

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GenePharma

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吉玛基因
www.genepharma.com

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1. Shen W , De Hoyos C L , Migawa M T , et al. Chemical modification of PS-ASO therapeutics reduces cellular protein-binding and improves the therapeutic index[J]. Nature biotechnology, 2019.
2. H , Lim S , Wong W F . ANTISENSE OLIGONUCLEOTIDES: FROM DESIGN TO THERAPEUTIC APPLICATION[J]. Clinical and Experimental Pharmacology and Physiology, 2010, 33(5-6):533-540.
3. Crooke S T. Molecular mechanisms of antisense oligonucleotides[J]. Nucleic acid therapeutics, 2017, 27(2): 70-77.
4. Cai Z, Cao C, Ji L, et al. RIC-seq for global in situ profiling of RNA–RNA spatial interactions[J]. Nature, 2020, 582(7812): 432-437
5. Li, M., Ding, X., Zhang, Y. et al. Antisense oligonucleotides targeting lncRNA AC104041.1 induces antitumor activity through Wnt2B/ β -catenin pathway in head and neck squamous cell carcinomas. Cell Death Dis 11, 672 (2020).

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